# INFRABED THERAPY FEATURES 6 JAN 22 Enlightened Therapies




# InfraBed Core Activation




Core Activation is the central therapeutic purpose\* for the InfraBed range of hyperthermal therapy products. Used for working both locally or generally\* to deeply penetrate for natural bio-activation with far infrared & subtle energies. \* "feeling and moving better"

# **Holistic Hyper-Thermia**



Holistic Hyperthermia is a most natural and popular holistic health practise. It includes Dry and Wet Saunas, Sun Bathing, Hot Springs, Hot Baths, Hot Lamps and Hyper Thermal Therapy InfraBeds. Controls are precise and easy. Add BioBlanket for a 'Cocoon'.

## **Amethyst Far Infrared Profile**



InfraBed's Amethyst grade A crystals, when heated, emit a subtle and deeply penetrating Far Infrared emission profile, around 9.2 numens. This is the 'sweet spot' considered ideal in bio-energetic thermal therapy. Infrared is natural, light-based, and anti-inflammatory in effect.

#### **IMulti Layers Construction**



InfraBed 15 layer construction is designed to protect from external emissions, incl EMF. Also to produce and reflect both ideal subtle and Far Infrared **Bio-therapeutic Emissions.** Several of your Infrabed layers, including a negative ions layer, are also in Energy BioBlankets

### **Big 3 ET InfraBed Warranties**




Enlightened Therapies takes its 'Big 3 Guarantee' : seriously. 1 is First Quality Manufacture 2 is Three Year Warranty and 3 is Complete Satisfaction so: ET guarantees you significant improvement of your issue or a full refund\* (within 40 days\*) Ph toll free ET 1300 305 107

## InfraBed Hyper Thermal Therapy




InfraBed Hyper-Thermal Therapy gives you influence on your Core Body Temperature with InfraBed's Time and Temperature controllers. These you can set and forget from 30 - 70 degrees C Temp and Timer from 3 to 8 hours, depending on types of ACHE.

#### InfraBed ACHE Settings



InfraBed has a range of setting options for your therapy needs. A cute use high-temp low-time. C hronic use mid-temp & time. Holistic low-temp & long-time. *E* nergetic room-temperature ie body temperature for any length of time, for holistic bioenergy protection and calming.

#### InfraBed 'Sun' Power



InfraBeds emit a natural 'Sun-like' Far Infrared, lightbased deep thermal therapy. Viruses and cancers do not like heat. Good core body temperature activates and supports organs, circulation, metabolism and immune system functionality.

## **Build Quality Certifications**



Your InfraBed range is Quality assured to CE standards and is certified by SAA approval for safety electrical standards in Australia and New Zealand. (240V, 288 Watt, Amps 1.22A) Guaranteed against any faults, in manufacture or materials with the ET full 3 Year Warranty.

## ET : Our Purpose



Colleen & Paul Keetley began ET as therapeutic device manufacturers with a purpose 'Helping Every Body Heal.' Non-invasive natural healing works this way so ET's goal is To Empower People to Safely Treat themselves and others.\* NB Consult your medical doctor

Empowerment